

Convoluting Ring Cushion

Unique central cavity for pressure free sitting

Designed to relieve pain in the tailbone and coccyx area



Convoluting surface suited to the heavier user



This cushion conforms to the contours of the body and is designed to relieve pain in the tailbone and coccyx area caused by pressure point discomfort. The convoluted surface improves ventilation to the skin's surface so ideal for people with increased skin sensitivity.

Made of one-piece, puncture-proof convoluted foam rubber. The unique design of the ring cushion means that vulnerable tailbone and coccyx region is effectively suspended by the central cavity for pressure free comfort when sitting. This helps minimise the risk of pressure injury when seated for long periods.

Ideal for post natal or post operative discomfort, haemorrhoids, hip bursitis, coccyx or prostrate inflammation; or other painful conditions.

Includes a lightweight stretch cover for added comfort, hygiene and ease of cleaning. The cover is easy to remove and can be machine washed.

SIZE	Large	Small
CODE	BL9535	BL9542
HEIGHT	90 mm	90 mm
WIDTH	370 mm	330 mm
DEPTH	450 mm	420 mm
CAPACITY	110 kg	110 kg
WEIGHT	160 g	120 g
WARRANTY	1 year	

