

Ideal for travel or temporary use, this secure raised toilet seat is easy to attach, remove and clean. Includes a flip up lid for improved appearance and increased hygiene.

The ergonomically designed seat has a tapered flange that sits directly on the porcelain bowl and features both a front and rear recess for increased hygiene access. With a smooth and gently contoured surface, the moulded one-piece seat is stain and odour resistant.

Three height options available to suit individual user needs or toilet heights.





CODE	BL1270	BL1263	BL1744
CAPACITY	130 kg	130 kg	130 kg
HEIGHT	Low 50 mm (2")	Medium 100 mm (4")	High 150 mm (6")
WEIGHT	1.25 kg	1.5 kg	1.8 kg
DIMENSIONS	407 mm deep x 370 mm wide		
	ARTG 176831		
CONFORMS TO	AS/NZS 3973:2009 , AS/NZS 4688.2:2000 & AS/NZS 4688.3:2000		
WARRANTY	1 year		

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH