

The Oxygen Bottle Carrier is designed to fit under the seat of any 'A' frame wheeled walker and will fit an oxygen cylinder up to 125 mm in diameter.

The holder is suspended by hooks that attach to the front and back cross piece of the walker and is set at an angle for optimal oxygen supply and to prevent the oxygen bottle from tipping out. It sits neatly underneath the seat of the walker.

Suitable for all BetterLiving wheeled walker models excluding the Tri-Wheel and Shopper.





CODE		BL1430	7
LENGTH		370 mm	_
HEIGHT		420 mm	_
MAX BOTTLE DIAMETER		125 mm	_
		ARTG 229461	
CONFORMS TO	370 mm 420 mm DIAMETER 125 mm ARTG 229461 AS EN 12182-200	AS EN 12182-2002 Technical Aids for Disabled Person	ns
WARRANTY		1 year	_
			_