

Ring Cushion



Unique central cavity for pressure free sitting

Designed to relieve pain in the tailbone and coccyx area



Designed to relieve pain in the tailbone and coccyx area, this ring cushion comfortably conforms to body contours to reduce pressure point discomfort.

Made of one-piece, puncture-proof latex foam rubber, the affected areas are suspended by the unique central cavity for pressure free comfort when sitting.

Includes a lightweight stretch cover for added comfort, hygiene and ease of cleaning. Machine washable at 30°C.

SIZE	18"	16"
CODE	BL9504	BL9511
HEIGHT	80 mm	80 mm
WIDTH	380 mm	320 mm
DEPTH	480 mm	410 mm
WEIGHT	0.4 kg	0.37 kg
CAPACITY	100 kg	100 kg
WARRANTY	1 year	



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.