## BetterLiving

Male Urinal Standard with Lid

## A lightweight urinal bottle with lid attached designed to help prevent spills.

Ideal for bed bound or immobilised males; and especially useful when there is little time between the urge to urinate and the ability to reach the toilet.

- o lid attached to the bottle; designed to help prevent spills
- lightweight translucent bottle, free of sharp edges inside and out

Graduated markings to 1000 ml

- graduated measurement markings up to 1000 ml, to accurately measure liquid quantity
- sturdy grip for easy hold and handling
- durable and easy to clean and disinfect (autoclave at 135°C or sterilise in boiling water or with steam)
- can be used when seated or standing, in several alternative positions

CONFORMS TO ARTG 202562 WARRANTY 1 year DUE TO HYGIENE REGULATIONS PRODUCT CANNOT BE RETURNED AFTER PURCHASE

103 g

BL1287

1000 ml

250 x 100 x 90 mm

High Density

Polyethylene

CODE

CAPACITY

DIMENSIONS

MATERIALS

WEIGHT

Sturdyhandle

Attachediloto prevent spills and reduce odour

