

Weighing just over 6 kilograms, the stylish Euro Lightweight Wheeled Walker is constructed from lightweight aluminium and designed for easy, every day use.

Height adjustable to best match the user, handle height is easily adjustable with an extension range of 100 mm to suit users of average height and weight. Moulded handgrips with soft touch brakes are ideal for people with limited strength.

Features a soft, wide nylon seat and stable backrest for maximum comfort when seated. Also includes a large, detachable, weatherproof fabric bag (with removable shoulder strap) that can hold up to 4.5 kg of shopping or other personal belongings; and a stick holder to ensure safe storage and transportation of walking stick.

This stylish European designer piece folds effortlessly with only one hand - simply pull the handle to collapse the walker for compact transportation or storage.













815 mm min 945 mm max



780 mm



615 mm



250 mm



450 mm

SEAT HEIGHT



540 mm

FOLDED WIDTH



270 mm



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL