

Softtouchbrakingmeet

Heightadiustablehandles

Lightweight alloy with tube in tube Lynweignen for increased strength

folddowndesign

With lower overall height range and seat height, the Petite Wheeled Walker is ideally suited for users 166cm (5'5") or less.

The alloy frame makes this wheeled walker especially lightweight, while the tube-in-tube construction significantly increases stability and durability.

Recommended for general indoor use, six inch castors ensure easy operation and manoeuvrability on smooth surfaces and through doorways.

The durable moulded handgrips with soft touch brakes are ideal for people with limited strength or mobility; maximising user comfort by minimising pressure on the hand. Handle height can be adjusted easily to best suit the user, with an extension range of 135 mm to suit users with a more petite stature.

The fully moulded PU seat is completely seamless, reducing the potential for skin tear or abrasion, plus foam backrest for added comfort when seated.

Cleverly concealed beneath the seat is a PVC pouch, offering additional security for small personal valuables.



Durable moulded orthopedic handgrip

Foambackrest

houlded seamless seat with Concealed PNC seat Pouch

Removable basket and bag

Idealfor people shorterthan

166 cm (5'

BetterLiving


With the added convenience of both basket and bag that can be used together (basket fits inside the bag to conceal contents) or individually. Both items are removable, for the security and transport of goods when leaving the wheeled walker unattended.

Easy closing mechanism and compact fold down design for convenient transport or storage.

COLOUR	Cherry Red
CODE	BL1126
CAPACITY	130 kg
WEIGHT	7.6 kg
CONFORMS TO	ARTG 176830
WARRANTY	1 year

HEIGHT










SEAT WIDTH







850 mm min 715 mm max

597 mm

580 mm

310 mm

360 mm

470 mm

210 mm

153 mm 6 inch

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL

