

Ideal for negotiating tight spaces, the Tri-Wheel Wheeled Walker is designed for enhanced manoeuvrability and smaller turning radius, enabling users to navigate obstacles and doorways otherwise impassable with a standard four wheeled model.

The lightweight frame features an easy fold mechanism for transport and storage. Height adjustable handles include light touch handbrakes and lockable park brakes for added security when stopped or for transfers.

The large 203 mm (8") solid wheels are puncture free and work well on all surfaces. The Tri-Wheel includes a large vinyl bag and easily removable front basket with tray.

Folded, the walker fits easily in any boot or even the backseat of most cars.







Tri-Wheel Wheeled Walker



COLOUR	Red
CODE	BL1356
CAPACITY	100 kg
TRAY AND BAG CAPACITY	5 kg
WEIGHT	7.9 kg
CONFORMS TO	ARTG 176830
WARRANTY	1 year

HEIGHT




970 mm min 870 mm max





660 mm



610 mm



190 mm



230 mm



750 mm





200 mm





203 mm 8 inch

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL