

With maximum overall width and weight tested to 220 kilograms, the Maxi Plus Wheeled Walker is the ideal choice for users of 160 kilograms or more. The sturdy steel frame with tube-in-tube construction increases product strength and durability.

Recommended for general indoor and outdoor use, eight inch castors provide increased mobility on uneven, outdoor or rough terrain surfaces.

The durable moulded handgrips with soft touch brakes are ideal for people with limited strength or mobility; maximising user comfort by minimising pressure on the hand. Handle height can be adjusted easily to best suit the user, with an extension range of 130 mm to suit users of very large or heavy stature.

With maximised width to accommodate the larger user, the moulded PVC seat is completely seamless, reducing the potential









for skin tear or abrasion, plus foam backrest for added comfort when seated.

Cleverly concealed beneath the seat is a PVC pouch, offering additional security for small personal valuables.

With the added convenience of both basket and bag that can be used together (basket fits inside the bag to conceal contents) or individually. Both items are removable, for the security and transport of goods when leaving the wheeled walker unattended.

Easy closing mechanism and compact fold down design for convenient transport or storage.

COLOUR	Cherry Red
CODE	BL8248
CAPACITY	220 kg
WEIGHT	12.4 kg
CONFORMS TO	ARTG 176830
WARRANTY	1 year


HEIGHT



920 mm min 790 mm max

LENGTH

690 mm



705 mm



335 mm





500 mm



565 mm



250 mm



8 inch

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL