

Shoprider COMO

ASK FOR THE BRAND SHOPRIDER



SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

SHOPRIDER[®]
www.shoprider.com.au *mobility*

COMO

The SHOPRIDER™ **COMO** portable powerchair is a well featured powerchair that offers the user the ability to easily fold and transport by using quick disconnect levers in seconds. The **COMO** features a comfortable mid back seat, adjustable height armrests, easy to use joy stick control and when combined with the powerful motors it will deliver outstanding safety & comfort.

SPECIFICATIONS

SPECIFICATIONS	COMO
Length	725mm
Width	545mm
Height	1180mm
Battery Capacity	12 Amp
#Maximum Range	up to 16km
Maximum Load Single User*	113Kgs
Charger on board	2Amp
Motor Size	0.50hp
Maximum Speed	5-6Kph
*Safe Climbing Angle	6 Degrees
*Max Climbing Angle	6 Degrees
Turning Radius	685mm
Wheel Size:	
- Front Caster	2 x 155mm
- Rear Caster	2 x 2 mm

FEATURES

- 12amp SLA Battery
- Mid Back Captain Seat
- Variable seat height
- Adjustable arm rests
- Controller interchangeable
- Free wheel device
- Easy Dismantle for Transport
- Total Weight 47.4kg
- Largest weight 18.4kg in batts front Chassis


Driver weight may exceed weight of unit.
Reduce speed when turning.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

COMO



SHOPRIDER[®]
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au