


# Dania Coloured Toilet Seats



Available in three colours

option for seat and lid, or seat only

A range of Dania plastic toilet seats available in a variety of colours, ideally suited to dementia patients. Fitted onto any standard toilet pan to replace the standard seat, these colourful seats help to provide a colour contrast with the white pan or wall for people with dementia or vision impairment.

Contrast is accepted as one of the most important aspects in a dementia care support environment, and these coloured toilet seats can help people to more easily locate the toilet.


These ergonomically designed toilet seats are compatible with most toilet pans and suited to most users. The flat-shaped seat facilitates easier sideways transfer for wheelchair users.

The contoured toilet lid is comfortable to sit on when the toilet is not in use, and acts as a backrest when upright.

Available in three colours, with option for seat and lid or seat only.



# Dania Coloured Toilet Seats



ITEM	Toilet Seat with Lid	Toilet Seat without Lid	Toilet Seat with Lid	Toilet Seat without Lid	Coloured Toilet Seat with Lid	Coloured Toilet Seat without Lid
CODE	EP-R37035B84	EP-R27035B84	EP-R37112B84	EP-R27112B84	EP-R37108B84	EP-R27108B84
COLOUR	Red	Red	Anthracite	Anthracite	Blue	Blue
DIMENSIONS	w 370 mm x L 450 mm					
MOUNTING POINTS	120 > 170 mm					
CAPACITY	250 kg					
WARRANTY	3 year					



Blue