

# Etac Fresh bath board

Combines comfort and safety.

## The original

The unique design of Etac Fresh has set the standard for bath boards. It's front edge, placed on the bath's outer side is wide enough to provide maximum support for transfers. The mid-section is narrowed off to allow easy access for intimate hygiene. Fresh is available in two lengths to ensure a perfect fit to various sizes of bath.

## Secure mounting

The edge stops are easy to mount and secure by using the included plastic coin. Four friction linings keep the bath board securely on the bathtub.



*Etac Fresh has a convenient place for soap and shampoo next to the handle.*

## Features



### Edge stop

The edge stops are easy to adjust and sit firmly to ensure a safe transfer.

## Accessory



### Support handle

81600040  
The handle ensures a comfortable and secure grip.

### Description

Fresh, 69 cm (27¼")	81600011
Fresh, 74 cm (29¼")	81600021
Fresh, 69 cm (27¼") with handle	81600014
Fresh, 74 cm (29¼") with handle	81600024
Fresh, 69 cm (27¼") with handle retail packaging	81600114

### Dimensions

Length: 69 or 74 cm (27¼ or 29¼")  
Width: 36–27–30 cm (14¼–10¾–11¾")  
Thickness: 4 cm (1½")  
Adjustability, edge stops:  
Fresh 69 cm (27¼"): max. 63.6 cm (25") min. 39 cm (15¼")  
Fresh 74 cm (29¼"): max. 68.6 cm (27") min. 43.2 cm (17")

### Material

Polypropylene  
Edge stops: Polypropylene with TPE surface

### Weight

Fresh 69 cm (27¼"): 2.1 kg (4.1 lbs)  
Fresh 74 cm (29¼"): 2.3 kg (5.1 lbs)



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.

