Etac Optima toilet arm support

Provides a firm and comfortable grip.

Etac Optima toilet arm support

Etac Optima arm support has a soft, rounded end which provides a comfortable grip. The supporting leg has an adjustable foot which ensures perfect floor contact. The large wall plate ensures a strong fix to the wall.



Etac Optima Lockable toilet arm support


Etac Optima Lockable toilet arm support, locks in either the folded up or down position. It provides stable support when standing. The unique locking mechanism is easy to operate with one hand. The arm support has a soft, rounded end which provides a comfortable and secure grip.



80303001

80303004

80303007



Description

Optima (one piece) Optima lockable (one piece) Optima lockable, extended (one piece)

Dimensions

Length arm supports: 71.5 cm (28¼") Length arm support, extended: 91.5 cm (36") Floor to arm supports: 80 cm (31½") Floor to arm supports, folded up: 147.5 cm (58") Wall bracket: 16.5x30.4 cm (6½x12")

Weight

2.4 kg (5.3 lbs) Extended: 3.34 kg (7.4 lbs)

Material

Wall bracket: Powder-painted aluminium Arm supports tubing and supporting leg: powder-painted steel tubing Arm supports cover: Polypropylene







Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85° C / 185° F if so required.

Accessories



Attachment for signal device 80303027 Fits Optima and Optima Lockable



Extension 10 cm (4"): 80303030 20 cm (7¾"): 80303031 30 cm (11¾"): 80303032 Fits Optima and Optima Lockable.



Toilet paper holder 83030001 Fits Optima and Optima Lockable.



Supporting leg 80303040 For relieving strain on walls. Fits Optima Lockable only.



Floor pillar 80303028 For floor-mounting, arm supports height: 80 cm (311/2") Fits Optima and Optima Lockable.

