

# Etac Supporter toilet arm supports

Stable toilet arm supports mounted directly on the toilet.



## Functional arm supports

Etac Supporter arm supports are foldable and have comfortable rounded hand grips. With its neat design, Supporter blends naturally into the bathroom.

## Easy to attach

Supporter is simple to mount without the need for tools. It fits both floor and wall mounted toilets. If space around the toilet is limited, Supporter can be used with only one arm support.

*Supporter includes seat and lid.*

## Accessories

See "Accessories for Hi-Loo and Supporter" on page 50



### Description

Supporter

80303112-2

### Dimensions

Total width: 60 cm (23½")

Width between arm supports: 49 cm (19¼")

Arm supports height above top of seating: 25 cm (9¾")

Arm supports length: 41.5 cm (16¼")

### Weight

5 kg (11 lbs)

### Material

Frame and arm supports: powder coated tubular steel

Arm supports cover: polypropylene

Lid and seat: Polypropylene

Mountings: aluminium



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.