Etac Swift Commode

A three in one product, that can be used as a commode, shower chair or toilet seat raiser.

Comfortable and functional

Etac Swift Commode is comfortable and discreet in the bedroom. The seat pad conceals a removeable pan with a secure lid that locks automatically when lifted. Swift Commode can be set with a slight anterior tilt to facilitate standing up. The soft seat and back support paddings come as standard.

3 in 1

Swift Commode can also be used as a shower chair or as a freestanding toilet seat raiser. The soft ferrules give an excellent grip on all surfaces.

Easy cleaning

Swift Commode is designed with easy cleaning in mind. It has smooth surfaces, easy detachable padding and is non corroding. Important features that simplify cleaning. Swift Commode comes flat packed and is easy to assemble and height adjust, without using tools.



With the padding, Swift commode has a discrete design and is comfortable to use when dressing and undressing.



Description

Swift commode

81702030

Dimensions

Total width: 56 cm (22") Seat width: 54 cm (211/4")

Seat height 42-57 cm (161/2-221/2")

Width between the arm supports: 45 cm (17¾")

Material

Seat, back and arm supports: polypropylene

Legs: aluminum and polyamide

Ferrules: TPE

Back and seat pads: Polythylene with nylon fabric

Weight

5.4 kg (12 lbs)











Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.



Integrated, self locking lid

The pan with a secure locking lid, fits smoothly and becomes one with the seat.

Accessory



Soap dish / Lid holder

80209266

Also handy for toilet paper Material: polypropylene

