

Etac Swift shower stool/chair

Swift is a prize-winning series of shower stools/chairs designed with functionality and simplicity in mind. The basic shower stool can easily be transformed to a complete shower chair, by adding arm- and back supports.

Stable on any floor

Etac Swift has a unique design that adapts to uneven floors thanks to the flexible seat. The soft, high quality ferrules adapt to the floor and ensure a good grip. These features contribute to the users safety.

Changes with the circumstances

Swift is adaptable to different needs. It is possible to start using the shower stool and add arm supports and back support if the circumstances change. Swift is easy to adjust to the required height with the telescopic legs. It can be set one step lower at the front for a slightly forward tilted seat to facilitate standing up. Soft pads which are specially suited for showering are available should extra comfort be needed.

Easy to handle

Swift is easy to handle thanks to its light weight. It comes flat-packed and is assembled without the need for any tools. Swift is also corrosion free and comes with a 5 year warranty!



Etac Swift

	Green	Grey	Blue
– shower stool	81701010	81701410	81701510
– shower chair with back and arm supports	81701030	81701430	81701530
– shower stool Low		81701450	
– shower chair Low w. back and arm supports		81701460	

Dimensions

Total height with back supports: 80–95 cm (31½–37½")
 Total height Swift Low with back supports: 72.5–80 cm (28½–31½")
 Total width, stool: 54 cm (21¼")
 Total width, with arm supports: 54.5 cm (21½")
 Total depth: 50 cm (19¾")
 Total depth Swift Low: 46.5 cm (18")
 Width between arm supports: 45 cm (17¾")
 Seat height: 42–57 cm (16½–22½")
 Seat height Swift Low: 34.5–42 cm (13½–16½")
 Seat area: 54x41cm (21¼ x 16¼")

Material

Polypropylene, aluminium

Weight

Swift shower stool: 3.1 kg (6.8 lbs)
 Swift Low shower stool: 2.9 kg (6.4 lbs)
 Swift shower chair: 4.7 kg (10.4 lbs)
 Swift Low shower chair: 4.5 kg (8.8 lbs)

Colours



red dot



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.

Features



Easy to fit
The push button lock function makes it easy to both fit and remove the legs.



Ferrules with a firm grip
The ferrules squeeze away the water and adapt perfectly to the floor.



Firm grip possibilities
Without arm supports, the openings enable a stable grip.



Safe surface
The surface structure ensures a safe, non slip seat, even when wet.



Etac Swift shower stool

Etac Swift shower chair with back and arm supports

Etac Swift shower stool Low with back and arm supports

Accessories



Back support
Green: 81701120
Grey: 84005069
Blue: 84005071



Back support pad
Green: 81701130
Grey: 84005072
The non-slip pad provides an instantly warm feeling and is water resistant. Material: polyethylene. Cleaning at max. 60°C (140°F)



Soap dish
Green: 81701140
Grey: 80209266
Blue: 81701340



Seat pad
Green: 81701132
Grey: 84005073
The non-slip pad provides an instantly warm feeling and is water resistant. Material: polyethylene. Cleaning at max. 60°C (140°F)



Arm supports pair
Green: 81701110
Grey: 84005074
Blue: 84005075

