

Etac Flex grab rail

Innovative design that creates possibilities.



Flexible modular grab rails

Etac Flex grab rail is the perfect choice for those who need a little extra support in the bathroom. Thanks to its flexible modular function, it is easy to adjust the handle to suit both you and your bathroom environment. The pure Scandinavian design blends in with other products and you can choose between two colours - white and grey.

Support for good and bad days

The material of the grab rail and its welldesigned shape and surface provide a stable grip even if your hands or the handles are wet. The triangular shape and the thumb support also make you feel safer.



Etac Flex fits perfectly in the bathroom environment.

Etac Flex screw mount		Weight
- Grab rail Flex 30 cm white	81707110	0.35 kg per unit (0.77 lbs)
- Grab rail Flex 60 cm white	81707120	
- Grab rail Flex 90 cm white	81707130	Maximum user weight
- Extension Flex 30 cm white	81706110	100 kg
- Grab rail Flex 30 cm grey	81707310	
- Grab rail Flex 60 cm grey	81707320	Material
- Grab rail Flex 90 cm grey	81707330	Polypropylene
- Extension Flex 30 cm grey	81706310	
		Colour
Etac Flex adhesive mount		White and Volcano grey
- Grab rail Flex 30 cm white	81707210	
- Grab rail Flex 60 cm white	81707220	
- Grab rail Flex 90 cm white	81707230	
- Extension Flex 30 cm white	81706210	
- Grab rail Flex 30 cm grey	81707410	3 min
- Grab rail Flex 60 cm grey	81707420	Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max 85°C / 185°F
- Grab rail Flex 90 cm grey	81707430	if so required.
- Extension Flex 30 cm grey	81706410	•

Note; To add one or more grab rail to a 30 cm, 60 cm or 90 cm combination, it is necessary to use the Etac Flex Extensions.

V

Features



Easy to adjust

The modular system can easily be assembled into a length and shape that suits you and your bathroom.



Secure grip

The well-designed shape and surface of the grab rail provide a stable grip even if your hands are wet.



Creating Possibilities

You can connect as many grab rails as you want and choose the angles that suit you. This creates possibilities.



Attach in two ways

Attach with screws or adhesive. The latter is perfect for those who cannot or do not want to drill in the bathroom.





Etac Flex offers a wide range of positioning options when connecting grab rails. Choose the angle between -90° to +90° that suits your need.

