

Encourages correct posture;  
provides excellent lumbar and lateral support

Available in  
Standard, High Back or with Seat



The Therapod back support system is versatile, durable and effective in relieving and preventing back pain in situations where alternative seating is not available. Designed to be completely adjustable and flexible enough to support every back, it is portable and lightweight and can be used in the car, home or at work.

It has a unique shell shaped structure with a series of adjustable straps that maintain full back support, which work by relieving lumbar disc pressure and reducing muscle activity in the back, providing pain relief for those who suffer from chronic back pain.

The shell design fits any seat and is finished with a cover of high grade urethane, upholstered in a soft and durable fabric.

The unique adjustment system allows the angle of the back rest to be set as required (an angle greater than 90° is recommended). Easy to install with four elasticised straps behind the zipper which are adjusted to vary the pressure at various points on your back.

Is fitted with a strap to attach to chair backs. Available in Standard or High Back, with optional coccyx and neck support.

Designed by an Australian Physiotherapist and Orthotist.

### Therapod® Coccyx Seat

The Therapod Coccyx Seat relieves sciatic and coccyx discomfort, helps stabilise the pelvis by tilting it forward and improves seat and back pressure distribution, particularly in car seats.

The seat can be used on its own or easily attaches to the Standard and High Back Therapod back support cushions.



therapod™

# Support Cushions



PORTABLE, LIGHTWEIGHT AND FULLY ADJUSTABLE TO SUIT THE INDIVIDUAL

ITEM	Standard	High Back	Seat	Neck Support
CODE	20BC01	20BC02	20BC00	20BC13
HEIGHT	600 mm	800 mm	50 mm	160 mm
WIDTH	500 mm	500 mm	500 mm	240 mm
DEPTH	100 mm	100 mm	500 mm	80 mm
WEIGHT	1kg	1kg	500 g	100 g
COLOUR	Black			
CONFORMS TO	ARTG175006			
WARRANTY	1 year			

