

Hematite Lightweight Transit Wheelchair

RG4600



<u>Features</u>

- Lightweight aluminium construction
- Super Lightweight 8.0kg Net Weight
- Powder coated folding frame
- Padded arm rests
- Fold down backrest
- 8" Solid rear tyres
- 8" Front rubber tyres
- Hard wearing nylon upholstery
- Swing away leg rests
- Non skid foot plates with heal strap
- Manual brake on rear wheels
- Seat belt fitted
- Rear storage pocket on backrest
- Available in marble finish Red

<u>Dimensions</u>

- Seat Width 18" (460mm)
- Seat Depth 400mm
- Seat Height 480mm
- Max User Capacity 110kg
- Wheelchair Weight 8.0kg

Options Elevating Legs Rest



E-mail: orders@shoprider.com.au



REDGUM Brand is associated with Shoprider Mobility www.redgumbrand.com.au or www.shoprider.com.au Ask for the Brand SHOPRIDER