

Sapphire Deluxe Lightweight Aluminium Wheelchair

RG18ALD / RG20ALD





Features

- Lightweight aluminium construction
- Powder coated folding frame
- Padded arm rests
- Lift up armrests
- 24" Solid tyres
- 7" Front rubber tyres
- Height adjustable seat
- Hard wearing nylon upholstery
- 50mm comfort cushion
- East to operate swing away leg rests
- Manual brake on rear wheels
- Easy remove rear wheels
- Available in Black or Red Sparkle finish

Dimensions

- Seat Width 18" (470mm) RG18ALD
- Seat Width 20" (510mm) RG20ALD
- Seat Depth 465mm
- Seat Height 380mm(+50mm cushion)
- Max User Capacity 115kg
- Wheelchair Weight 18" 13.5kg
- Wheelchair Weight 20" 14.0kg

Options Available

• Anti tip wheels

Phone: 08 9248 4180 Fax: 08 9248 9452

E-mail: orders@shoprider.com.au



REDGUM Brand is associated with Shoprider Mobility www.redgumbrand.com.au or www.shoprider.com.au Ask for the Brand SHOPRIDER