

Invacare® | **Colibri™**

**NEW**  
with  
18Amp Battery



**LiteLock™**  
system  
(patent pending)



## The colourful route to effortless, safe drive

The new *Invacare Colibri* is a stylish, simple and colourful micro lightweight scooter, designed for those who enjoy a totally independent lifestyle.

Thanks to the unique *Invacare LiteLock™* system the scooter can easily be taken apart without the need for tools and fits neatly into any car boot. Plus, with its small footprint the *Invacare Colibri* offers greater manoeuvrability in tight, or limited spaces. That's why the *Invacare Colibri* goes anywhere your fancy takes you!

Additionally, you can choose from a range of trendy colour panels to match your style, swapping them according to your needs.

Be ready to roll comfortably and securely seated in your new stylish companion!



**Yes, you can.®**



## Features and Options



Choose from five great colours

### LiteLock system - easy to take apart and reassemble

The ingenious and unique *Invacare LiteLock™* system automatically releases once the battery box has been removed. It takes less than a minute to take the scooter apart – without tools. Its front section can be simply lifted and placed in the boot of a car, along with the other sections. Later, it's just as easy to reassemble and the scooter is secure.



### Effortless dismantling

Don't think - do it! Less than a minute is all that's needed to disassemble the **Colibri** for transportation.



### Increased comfort

The soft, padded seat provides comfort during driving. For increased comfort there is also the optional 'easy fit' seat suspension kit.



### Ergonomic Handlebar and Plastic Basket

For greater comfort and prevents fatigue in hands and wrists.

## Technical data

Invacare <b>Colibri</b>	465 mm	400 mm	49cm	385 mm	505 mm
Invacare <b>Colibri</b>	1010 mm	49 kg (18 Ah batteries)	136 kg	2 X 18 Ah	2200 mm
Invacare <b>Colibri</b>	45 mm	6° / 10.5%	200 W	Up to 16 km	8 km/h



Invacare Australia Pty Ltd  
1 Lenton Place (PO Box 5002)  
North Rocks NSW 2151, Australia  
Freephone: 1800 460 460  
Fax: 1800 814 367  
Email: sales@invacare.com.au  
Web: www.invacare.com.au

Invacare New Zealand  
4 Westfield Place (PO Box 62-124)  
Mt Wellington Auckland, New Zealand  
Freephone: 0800 INVACARE (0800 468 222)  
Freefax: 0800 807 788  
Email: sales@invacare.co.nz  
Web: www.invacare.co.nz

## Shroud colours

Scooter Colibri white 4w	Scooter Colibri jasper red 4w	Scooter Colibri jade green 4w	Scooter Colibri sapphire blue 4w	Scooter Colibri stone orange 4w
pkg-1577168	pkg-1577169	pkg-1577170	pkg-1577171	pkg-1577172

1/ Compliant with the directive on the Restriction of the use of certain Hazardous Substances in electrical and electronic equipment - RoHS (2011/65/EC)



Yes, you can.®

© 2014 Invacare International Sàrl  
All rights reserved.

PB6835.05.14