

P321

## Ezi-go

### FEATURES

- **115kg Seating Capacity.**
- Dismantles in seconds for easy transportation and storage
- Puncture proof tyres front & rear.
- Quick release seat and armrests
- Armrests adjust for height and width
- Suitable for indoors and shopping centres
- Latest Technology "PG Pilot" Joystick Controller.
- Joystick mounting for left and right hand operation
- Quick release trigger separates the base sections.
- Thickly padded seat gives greater comfort to regular users.
- Intelligent, regenerative electro magnetic brakes activate automatically when Joy Stick is released



### Dimensions

Maximum Speed	6.5kph
Maximum Range	5-8km
Motors	2 x 120 watt
Batteries	2 x 12 volt sealed gel type
Charger	1.5 amp "Switch Mode"
Gradient	8 deg
Front Caster Wheels	150x50 Foam Filled
Rear Drive Wheels	200x75 Foam Filled
Length (footplate up / down)	70 - 87 cm
Width (armrests)	56 cm
Height (with backrest in normal position)	79 - 84 cm
Seat Width	46 cm
Seat Height From Ground	49 - 54 cm
Seat Height From Foot Plate	39 - 44 cm
Weight of Seat inc armrests and Joystick	12 kg
Total Weight including Batteries	41 kg