

Scooter Accessories

Scooter Xtras has a range of accessories to assist in the daily use of mobility equipment. These accessories have been designed to deliver value for money and functionality for the user of the mobility Scooter.

Bag Carrier - SE-BC

Colour – Black.

Frame – Powder Coated Steel.

Bag – Waterproof Vinyl with hard board base.

Capacity – 12kg (10kg when fitted with Sun Canopy).

Features

Integrated Walking Stick / Crutch Carrier.

Rear mounting tube can be cut to suit most scooter models.

Sun Canopy may be added as option Ready to assemble.

*Walking stick not included



Deluxe Sun Canopy - SE-SD

Colour – Red, Black or Steel Blue

Frame – Powder Coated Steel.

Fabric – Poly Cotton .

Canopy Dimensions – 950mm long by 650mm wide.

Features

Easy 200mm Height Adjustability Folding Side Flaps for Added Protection from the Sun





Scooter Xtras is a Shoprider Australia Pty Ltd Brand www.shoprider.com.au



Sun Canopy - SE-S New May 2014

Colour – Red or Black.

Frame – Powder Coated Steel.

Fabric – Poly Cotton .

Canopy Dimensions – 770mm long

by 560mm wide.

Features

Sleek new design

Larger mid size

Height adjustable 120mm

Ready to assemble.

Folding Side Flaps for Added



Wheelie Walker Carrier -SE-WWB

Colour - Black

Frame – Powder Coated Steel.

Features

Rear mounting tube can be cut to suit

most scooter models.

Suits most seat walkers

Integrated Walking Stick / Crutch

Carrier.

Ready to assemble.

*Walker not included





Scooter Xtras is a Shoprider Australia Pty Ltd Brand www.shoprider.com.au



Walker Carrier – SE-WFC

Colour – Black

Frame – Powder Coated Steel.

Features

Rear mounting tube to suit most scooter models.

Suits most walking frames Ready to assemble.

Now in packaged in a carton for easy transport.

Can be Air Freighted



Walking Stick Holder – SE-WSH

Colour – Black

Frame – Powder Coated Steel.

Features

Rear mounting tube to suit most scooter models.


Can be mounted on side of the seat on some Shoprider models

Ready to assemble.

*Walking stick not included

Now in packaged in a carton for easy transport.

Can be Air Freighted



Scooter Xtras is a Shoprider Australia Pty Ltd Brand www.shoprider.com.au



Oxygen Bottle Carrier - TE-O

Colour – Black

Frame – Powder Coated Steel.

Features

Rear mounting tube to suit most scooter models.

Can be mounted on side of the seat on some Shoprider models Integrated Walking Stick / Crutch Carrier. *Walking stick not included

Ready to assemble.



Flag Holder - SE-FH

Colour – Black

Frame - Powder Coated Steel.

Features

Rear mounting tube to suit most scooter models.

Can be Mounted on side of the seat on some Shoprider models

Ready to assemble.

* Flag not included

Can be Air Freighted



All of the above accessories come complete with all the necessary fitting required to fit the accessory securely to the under seat mounting point. The under seat mounting point needs to be 25mm x 25mm to accommodate these accessories.

For your nearest Scooter Xtras distributor Call 1800 672 633

Scooter Xtras is a Shoprider Australia Pty Ltd Brand www.shoprider.com.au

Your Nearest Authorised Distributor

Correct as 1/5/2017 and subject to change without notice