

## **Pearl Deluxe Aluminium Wheelchair**

Part Number RG18AAD / RG20AAD / RG18AADRHB / RG20AADRHB



## **Features**

- Aluminium anti corrosive folding frame
- Padded adjustable height arm rests
- Adjustable height rear handles
- Adjustable comfort back rest (Velcro)
- 24" Solid tyres
- 7" Front rubber tyres
- Adjustable castor angle
- 4 stage seat adjustment
- Hard wearing nylon upholstery
- Height & angle adjustable foot plates
- Easy to operate lever swing away leg rests
- Manual brake on rear wheels
- Quick disconnect rear wheels
- Anti tip wheels fitted as standard
- Available in Pearlescent White

## **Dimensions**

- Seat Width 18" (470mm) 20" (510mm)
- Seat Depth 440mm
- Seat Height 410mm
- Max User Capacity 18" 120kg
- Max User Capacity 20" 130kg
- Wheelchair Weight 18" 17.0kg
- Wheelchair Weight 20" 18.5kg

<u>Rear Hand Brake Model</u> <u>now available</u>

Phone: 08 9248 4180 Fax: 08 9248 9452 E-mail: orders@redgumbrand.com.au



**REDGUM Brand is a 100% Australian Owned & Operated**