

## **Updated Model Deluxe Folding Transporter Wheelchair**





The REDGUM Brand is owned and operated by Shoprider Australia Pty Ltd. <a href="https://www.shoprider.com.au">www.shoprider.com.au</a>

## **Features**

- Lightweight aluminium construction powder coated folding frame
- Fold down back rest
- Padded arm rests
- Lift up armrests
- Removable armrests
- 12" Rear with durable rubber
- 8" Front rubber tyres
- Hard wearing nylon upholstery
- Swing away leg rests
- Attendant rear control handbrakes
- Manual brake on the rear wheels
- Seat belt fitted
- Rear storage pocket on backrest
- Available in Red & Blue

## **Dimensions**

- Seat Widths 17"(420mm)19"(470mm) or 22" (550mm)
- Seat Depth 390mm
- Seat Height 460mm
- Max user capacity RG17HBD & RG19HBD -110KG
- Max user capacity RG22HBD -125kg
- Wheelchair weight RG17HBD & RG19HBD -13.5kg
- Wheelchair weight RG22HBD -16.5kg

Your Nearest Authorised Distributor