

Scooter Accessories

Scooter Xtras has a range of accessories to assist in the daily use of mobility equipment. These accessories have been designed to deliver value for money and functionality for the user of the mobility Scooter.

Bag Carrier - SE-BC Colour – Black. Frame – Powder Coated Steel. Bag – Waterproof Vinyl with hard board base. Capacity – 12kg (10kg when fitted with Sun Canopy). <u>Features</u> Integrated Walking Stick / Crutch Carrier. Rear mounting tube can be cut to suit most scooter models. Sun Canopy may be added as option Ready to assemble. *Walking stick not included



Deluxe Sun Canopy - SE-SD Colour – Red, Black or Steel Blue Frame – Powder Coated Steel. Fabric – Poly Cotton . Canopy Dimensions – 950mm long by 650mm wide. Features Easy 200mm Height Adjustability

Folding Side Flaps for Added Protection from the Sun









Sun Canopy - SE-S Colour – Red or Blue . Frame – Powder Coated Steel. Fabric – Poly Cotton . Canopy Dimensions – 620mm long by 440mm wide. <u>Features</u> Foldable to the rear . Upright tube can be cut to suit the correct height.. Ready to assemble. *Walking stick not included



Wheelie Walker Carrier -SE-WWB Colour – Black Frame – Powder Coated Steel. <u>Features</u> Rear mounting tube can be cut to suit most scooter models. Suits most seat walkers Integrated Walking Stick / Crutch Carrier. Ready to assemble. *Walker not included





Scooter Xtras is a Shoprider Australia Pty Ltd Brand www.shoprider.com.au



Walker Carrier – SE-WFC Colour – Black Frame – Powder Coated Steel. <u>Features</u> Rear mounting tube to suit most scooter models. Suits most walking frames Ready to assemble. Now in packaged in a carton for easy transport. Can be Air Freighted



Walking Stick Holder – SE-WSH Colour – Black Frame – Powder Coated Steel. <u>Features</u> Rear mounting tube to suit most scooter models. Can be mounted on side of the seat on some Shoprider models Ready to assemble. *Walking stick not included Now in packaged in a carton for easy transport. Can be Air Freighted



Scooter Xtras is a Shoprider Australia Pty Ltd Brand www.shoprider.com.au



Oxygen Bottle Carrier - TE-O Colour – Black Frame – Powder Coated Steel. <u>Features</u> Rear mounting tube to suit most scooter models. Can be mounted on side of the seat on some Shoprider models Integrated Walking Stick / Crutch Carrier. *Walking stick not included Ready to assemble.



Flag Holder - SE-FH Colour – Black Frame – Powder Coated Steel. <u>Features</u> Rear mounting tube to suit most scooter models. Can be Mounted on side of the seat on some Shoprider models

Ready to assemble. * Flag not included Can be Air Freighted



All of the above accessories come complete with all the necessary fitting required to fit the accessory securely to the under seat mounting point. The under seat mounting point needs to be 25mm x 25mm to accommodate these accessories.

For your nearest Scooter Xtras distributor Call 1800 672 633

Scooter Xtras is a Shoprider Australia Pty Ltd Brand www.shoprider.com.au Your Nearest Authorised Distributor